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DNA Gel Extraction Kit Protocol

1. Excise less than 500 mg gel slice containing the DNA fragment with a clean scalpel or razor blade under ultraviolet lamp and transfer to a 1.5 ml microcentrifuge tube (not provided).

- * Minimize the size of the gel slice by removing extra agarose to reduce dissolve time.
- * The gel weight should not exceed 500 mg, otherwise it will block the spin column or affect the DNA recovery efficiency.

2. Add 500 µl Buffer G, incubate at 50°C until gel slice is completely dissolved(about 10 minutes).

- * The dissolution time will increase as the volume of the gel increases. For >2% agarose gels, more dissolution time is required.
- * Mix the tube by inversion every 2-3 minutes to help dissolve gel.
- * Check that the color of the mixture is yellow (similar to Buffer G without dissolved agarose). If the color of the mixture is violet, add 10 µl of 3 M sodium acetate, pH 5.0, and mix. The color of the mixture will turn yellow.
- * If more than 300 mg was used, it will be difficult to judge whether the gel is completely dissolved. Use a pipette to aspirate the liquid carefully and observe if there is any residual gel.

3. Add 200 µl isopropanol and mix thoroughly.

- * If the DNA fragment is between 500 bp and 4 kb, this step can be omitted.

4. Transfer the sample into the spin column (The spin column is placed in a 2 ml collection tube), Centrifuge at 12000 rpm for 30 s. Discard the filtrate and place the spin column back into the collection tube.

- * For sample volumes of more than 800 µl, it should be centrifuged through the column twice.
- * If the gel weight exceeds 300 mg or the concentration is $\geq 2\%$, maybe it needs to be centrifuged at full speed for 1 min again.

5. Add 500 µl Buffer WS to the spin column, closed the lid and centrifuge at 12000 rpm for 30 s. Discard the filtrate and place the spin column back into the collection tube.

- * This step is to remove the residual trace agarose molecules. If the recovered DNA is not for sequencing, in vitro transcription or microinjection experiments, this step can be omitted.

6. Add 700 µl Buffer WG, closed the lid and centrifuge at 12000 rpm for 30 s. Discard the filtrate and place the spin column back into the collection tube.

- * If the recovered DNA is used for salt sensitive experiments, such as blunt-end ligation experiment or sequencing, it is recommended to add Buffer WG and stand at room temperature for 2-5 minutes before centrifugation.
- * Ensure ethanol has been added into Buffer WG.

7. Repeat step 6 once.

8. Centrifuge at 14000 rpm for 1 min.

- * If the full speed of the centrifuge could not up to 14000 rpm, centrifuge at full speed for 2 min.
- * Do not pass over this step, or the residual ethanol in the eluted DNA will affects the final applications.

9. Discard the collection tube, place the spin column into a cleaning 1.5 ml centrifuge tube (provided in the kit), add 25-30 µl Buffer TE to the center of the column membrane, close the lid, stand at room temperature for 1 min, and centrifuge at 12000 rpm for 30 s.

- * The deionized water can also be used to elute DNA, but the pH value should be within 7.0 - 8.5, otherwise it may affect the elution efficiency.

10. Discard the column, the eluted DNA can be used for molecular biology experiment immediately, or store the DNA at -20°C.

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DNA Purification Kit Protocol

1. Add 5 times the volume of Buffer P to 1 volume of the PCR sample or the DNA solution that needs to be cleaned. Do not discard the tip, pipette several times to mix well, transfer the mixture to the spin column (the spin column is placed in a 2 ml collection tube), closed the lid.

* For example, add 500 µl Buffer P for every 100 µl PCR sample.

* If the PCR sample contains mineral oil, it does not need to be removed or included in the sample volume.

* DNA solutions for cleaning include enzymatic reaction solutions (such as enzyme digestion reactions, ligation reactions, etc.), RNase-treated DNA solutions (which can remove degraded RNA) and DNA containing impurities obtained by phenol/chloroform extraction.

* The dye added in Buffer P can indicate the pH. If the mixture turn to violet after adding Buffer P to the sample, add 10 µl of 3 M sodium acetate (pH 5.0) and mix, the color of the mixture will turn yellow, otherwise DNA adsorption will be inefficient.

2. Centrifuge at 12000 rpm for 30 s, discard the filtrate and place the spin column back into the collection tube.

3. Add 700 µl Buffer WB to the spin column, closed the lid, centrifuge at 12000 rpm for 30 s. Discard the filtrate and place the spin column back into the collection tube.

* Ensure ethanol has been added into Buffer WB.

4. Centrifuge at 14000 rpm for 1 min.

* If the full speed of the centrifuge could not up to 14000 rpm, centrifuge at full speed for 2 min.

* Do not pass over this step, or the residual ethanol in the eluted DNA will affects the final applications.

5. Discard the collection tube, place the spin column into a cleaning 1.5 ml centrifuge tube, add 30-50 µl Buffer TE to the center of the column membrane, close the lid, stand at room temperature for 1 min, and centrifuge at 12000 rpm for 30 s.

* The deionized water can also be used to elute DNA, but the pH value should be within 7.0 - 8.5, otherwise it may affect the elution efficiency.

6. Discard the column, the eluted DNA can be used for molecular biology experiment immediately, or store the DNA at -20°C.

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1. Huang Y, Feng H, Lu H, et al. Novel 16S rDNA primers revealed the diversity and habitats-related community structure of sphingomonads in 10 different niches[J]. Antonie van Leeuwenhoek, 2017, 110(7): 877-889.
2. Chen J, He Y, Wu Y, et al. Single Ethanol withdrawal regulates extrasynaptic δ -GABAA receptors via PKC δ activation[J]. Frontiers in molecular neuroscience, 2018, 11: 141.
3. Guo L, Song Y, Yuan Y, et al. Identification of nucleic acid aptamers against lactate dehydrogenase via SELEX and high-throughput sequencing[J]. Analytical and Bioanalytical Chemistry, 2021, 413(17): 4427-4439.
4. Wang P, Jia X, Xiao X, et al. An early diagnostic clue for COL18A1-and LAMA1-associated diseases: high myopia with alopecia areata in the cranial midline[J]. Frontiers in Cell and Developmental Biology, 2021, 9: 644947.
5. Yao Y, Deng S, Zhu F. Prenatal detection of novel compound heterozygous splice site variants of the KIAA0825 gene in a fetus with postaxial polydactyly type A[J]. Genes, 2022, 13(7): 1230.
6. Fu Y, Chen Q, Xiong M, et al. Clinical performance of nanopore targeted sequencing for diagnosing infectious diseases[J]. Microbiology Spectrum, 2022, 10(2): e00270-22.

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1. Zhang X H, Liu Z Q, Xue Y P, et al. Activity improvement of a regioselective nitrilase from *Acidovorax facilis* and its application in the production of 1-(cyanocyclohexyl) acetic acid[J]. Process Biochemistry, 2014, 49(12): 2141-2148.
2. Liu Z Q, Zhang X H, Xue Y P, et al. Improvement of *Alcaligenes faecalis* nitrilase by gene site saturation mutagenesis and its application in stereospecific biosynthesis of (R)-(-)-mandelic acid[J]. Journal of agricultural and food chemistry, 2014, 62(20): 4685-4694.
3. Huang X, Shi T, Mo K, et al. Monoclonal antibody against premembrane viral protein of avian tembusu virus[J]. Monoclonal Antibodies in Immunodiagnosis and Immunotherapy, 2017, 36(2): 57-61.
4. Huang X, Qiu H, Peng X, et al. Molecular analysis and serological survey of Tembusu virus infection in Zhejiang, China, 2010-2016[J]. Archives of virology, 2018, 163(12): 3225-3234.

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